



RBP Chemical Technology, Inc.
150 S. 118th Street • P.O. Box 14069
Milwaukee, WI 53214-0069
www.rbpchemical.com
Corporate: 414-258-0911
Fax: 414-258-7908
RBP East: 508-865-0800

Product Data Sheet

Date: 08/10/11
Supersedes: New
PRODUCT #: J0152

ANTISTAT SOLUTION

High-Performance Antistatic Concentrate

Description:

Antistat Solution is an additive to eliminate static build up in the paper and enhance the folding and stacking of the printed copies. This solution is designed to be added to the silicone solution to help this process. **Antistat Solution** contains additives and biocides for improving the silicone's performance.

Directions for Use:

Mix 4 to 8 ounces per gallon of **Antistat Solution** to each gallon of water. For silicone solutions, add directly to the mixed silicone utilizing a dosing pump or by hand at press side.

Specifications:

Appearance:	Opaque blue liquid
Density:	8.3 lbs/gal
VOC Content: (EPA Method 24)	None
Flash Point:	> 212 °F

This product should be used only for its intended purpose. The information stated above is based on our laboratory tests and experience, and is accurate to the best of our knowledge. Since actual use is beyond our control, the recommendations or suggestions are made without warranty, expressed or implied.