



RBP Chemical Technology, Inc.
150 S. 118th Street • P.O. Box 14069
Milwaukee, WI 53214-0069
www.rbpchemical.com
Corporate: 414-258-0911
Fax: 414-258-7908
RBP East: 508-865-0800

Product Data Sheet

Date: 08/01/11
Supersedes: New
PRODUCT #: J5203

AQUA SERIES 8500

Two Step Mild Acid Sheetfed Fountain Concentrate

Description:

A two step mild acid fountain concentrate designed for use on sheetfed presses. **AQUA SERIES 8500** is formulated with unique desensitizers and wetting agents to provide superior plate protection with excellent ink and water balance. It effectively addresses piling issues related to newer paper stocks and is compatible with all types of plates. **AQUA SERIES 8500** is a low dosage product.

Directions for Use:

Mix 2 oz/gallon of **AQUA SERIES 8500** to each gallon of water. Increase dosage if plate sensitivity is occurring and decrease dosage if ink emulsification is a problem. Each oz/gallon of **AQUA SERIES 8500** will increase conductivity by 900 micromhos/cm. Use **AQUA SERIES 8500** with an RBP alcohol replacer such as **AQUANOL[®] 600**, to eliminate the use of alcohol.

Specifications:

Appearance:	Clear blue green liquid
Density:	9.1 lbs/gal
Conductivity:	900 micromhos/cm per ounce
pH, diluted:	4.5-5.5
VOC Content:	0.82 lbs/gal, 98 g/L, 9.0%

This product should be used only for its intended purpose. The information stated above is based on our laboratory tests and experience, and is accurate to the best of our knowledge. Since actual use is beyond our control, the recommendations or suggestions are made without warranty, expressed or implied.

